



Purpose

In this experiment, we applied previous knowledge in reactions in order to identify two cations in an unknown solution.

Background

We gathered information about 9 different cations and their observable reactions to specific reagents.

1 of 3





Background Continued



2 of 3

Starts: clear and colorless Ends: cloudy white, white precipitate

Al(NO₃)₃ + NaOH

Starts: clear and colorless Ends: cloudy white, white precipitate, jelly-like

Starts: clear and yellow Ends: blood red, no precipitate

$$Mn(NO_3)_2 + NaBiO_3$$

Starts: Clear and colorless Ends: clear, forms brown precipitate

$$Ba(NO_3)_2 + Na_2SO_4$$

Starts: Clear and colorless Ends: pearl white, white precipitate



Background Continued



3 of 3

 $Ca(NO_3)_2 + (NH_4)_2C_2O_4$

Starts: Clear and colorless Ends: pearl white, white precipitate

Ni(NO₃)₂ + dimethyl glyoxime

Starts: Clear and light green Ends: pink/red, formed precipitate $Cu(NO_3)_2 + KI$

Starts: Clear and light blue Ends: yellow-brown, sludge-like

Co(NO₃)₂ + 1-nitroso-2-naphthol

Starts: Clear and pinkish red Ends: dark red, formed precipitate

Experimental Procedure

Experimental Procedure

Experimental Procedure

- Our solution began as a pinkish-red
- <u>Rxn 1:</u> HNO₃ and NH₃ were mixed with our solution, it turned beer like in color and did not form any precipitate. No aluminum, iron, or manganese is present. We can now focus on the right side of the flow chart.
- Rxn 2: (NH₄)₂CO₃ is mixed solution from rxn 1, becomes a grapefruit like color of pinkish red. A precipitate forms, which we centrifuge and decant. The liquid is a bright purple.
- Rxn 3: Our decanted liquid from rxn 2 is mixed with Na₂HPO₄, becoming a deeper purple, no precipitate is formed meaning it does not have magnesium.

- Rxn 9: HCl is mixed with solid
 b. It turns a very light pink
 color, meaning either Ba²⁺ or
 Ca²⁺ is in the substance.
- Rxn 10: Na₂SO₄ is mixed with the substance from rxn 9, it remains a clear color and no precipitate forms. This means that the substance contains Ca²⁺ and no Ba²⁺.
- Rxn 11: (NH₄)₂C₂O₄ is mixed with the substance from rxn 10. The solution forms a white precipitate, meaning our first cation is CaC₂O₄.

- Step 12: The liquid from rxn 3 is evaporated until dry, leaving a solid behind. The powder left is a purple color.
- Rxn 13: H₃PO₄ is mixed with the powder. The formed liquid is a light see through pink once the powder fully dissolves.
- <u>Rxn 14:</u> Na₂HPO₄ is mixed with the solution from rxn 13. It remains a see through pink with no precipitate.
- Rxn 15: KI is mixed with solution from rxn
 14. No change occurs, meaning it is not CuI.
- <u>Rxn 16</u>: C₄H₈N₂O₂ is mixed with solution from rxn 14. No change occurs, meaning it is not NiC₈H₁₄N₄O₄
- Rxn 17. C₁₀H₆(NO)OH is mixed with solution from rxn 14. It turns a deep blood like red, becomes sticky in texture, and forms precipitate. Our second cation is NaCO(NN)₃.



Summary & Conclusion

We utilized previous knowledge about the appearance and consistency of nine different cations when reacted with specific reagents to determine the identify of two unknown cations in an unknown solution. After following the flowchart and conducting a series of tests, we determined the possible presence of Ca²⁺ and CO²⁺ in our unknown solution.